



- Notes**
1. 标#号为重点管控尺寸
 2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
 3. 产品要求洁净免洗 生产场所 10万级净化车间
 4. 未标注尺寸参照3D档

STATE: Release	Select		V		UNITS	MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA								
		± A	± B	± C	MAT'L			PART NO. (INTENDED USE)	TITLE:							
	0-6	0.02	0.05	0.10	导电PP	FINISH	APPD:	罗氏801吸头								
	>6-30	0.05	0.10	0.15		O'TY	CHKD:	DWG NO.:								
	>30-60	0.10	0.15	0.20		DR:	<table border="1" style="display: inline-table; border-collapse: collapse;"> <tr> <td style="text-align:center;">SCALE</td> <td style="text-align:center;">SHEET</td> <td style="text-align:center;">REV.</td> </tr> <tr> <td style="text-align:center;">1:1</td> <td style="text-align:center;">1/1</td> <td style="text-align:center;">A</td> </tr> </table>				SCALE	SHEET	REV.	1:1	1/1	A
	SCALE	SHEET	REV.													
	1:1	1/1	A													
>60-120	0.15	0.20	0.25													
>120-180	0.20	0.25	0.30													
>180-300	0.25	0.30	0.35													