



A

A




B

B



C

C




细节 A  
比例 3.000

查看细节 A

D

D



## Notes

1. 标Φ为直径管径尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净无脱 生产场所 10万级净化车间
4. 表格尺寸参照3D图

STATE: Release

Select	v			UNITS	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA	
	± A	± B	± C	MM			
	0-6	0.02	0.05	0.10			MAT'L
	>6-30	0.05	0.10	0.15			PP
	>30-60	0.10	0.15	0.20			FINISH
	>60-120	0.15	0.20	0.25			APPD:
>120-180	0.20	0.25	0.30	O'TY	CHKD:	DWG NO.:	
>180-300	0.25	0.30	0.35	DR:	SCALE	SHEET	REV.
					1:1	1/1	A

E

A4