


REV.		DR:	DATE	NOTE	CRAFTWORK
------	--	-----	------	------	-----------

**Notes**

1. 标\*号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和磕伤
3. 产品要求洁净免洗 生产场所 10万级净化车间
4. 未标注尺寸参照3D档

STATE: Release

Select		V		UNITS	MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA								
	± A	± B	± C	MAT'L		PART NO. (INTENDED USE)	TITLE: PCR八连管盖								
0-6	0.02	0.05	0.10	PP		APPD:	DWG NO.:								
>6-30	0.05	0.10	0.15		FINISH		CHKD:	 <table border="1"> <tr> <td>SCALE</td> <td>SHEET</td> <td>REV.</td> </tr> <tr> <td>1:1</td> <td>1/1</td> <td>A</td> </tr> </table>			SCALE	SHEET	REV.	1:1	1/1
SCALE	SHEET	REV.													
1:1	1/1	A													
>30-60	0.10	0.15	0.20			DR:									
>60-120	0.15	0.20	0.25												
>120-180	0.20	0.25	0.30	O'TY											
>180-300	0.25	0.30	0.35												

E A4