


1	2	3	4
REV.	DR:	DATE	CRAFTWORK


A
A
B

A
B




C

C



D




Notes

1. 格号为重要管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净免洗 生产场所 10万级净化车间
4. 表格尺寸参照3D图

STATE: Release

Select			V	UNITS	MM	NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA
	± A	± B	± C	MAT'L			
0-6	0.02	0.05	0.10	PP	PART NO. (INTENDED USE)	TITLE:	100ul 双色全裙边
>6-30	0.05	0.10	0.15				
>30-60	0.10	0.15	0.20	FINISH	APPD:	DWG NO.:	
>60-120	0.15	0.20	0.25	Q'TY	CHKD:		
>120-180	0.20	0.25	0.30		DR:		
>180-300	0.25	0.30	0.35	1			



SCALE	SHEET	REV.
1:1	1/1	A

E

A4