


Notes

1. 标*号为重点管控尺寸
2. 产品表面要求光滑 允许不得有缺料、缩水、毛边、脏污和划痕
3. 产品要求洁净无油 生产场所 10万级净化车间
4. 表格尺寸参照3D档

STATE: Release

Select	UNITS			NAME (INTENDED USE)	浙江友莱生物科技有限公司 YOU LAI BIOTECHNOLOGY CO.,LTD. ZHE JIANG,CHINA				
	MM								
	± A	± B	± C		PART NO. (INTENDED USE)	TITLE: T300吸头			
	0-6	0.02	0.05						0.10
	>6-30	0.05	0.10		0.15	APPD:	DWG NO.:		
	>30-60	0.10	0.15		0.20				
>60-120	0.15	0.20	0.25	CHKD:	SCALE SHEET REV.				
>120-180	0.20	0.25	0.30						
>180-300	0.25	0.30	0.35	DR:	1:1	1/1	A		



E

A4