

1

2

3

4



REV.

DR:

DATE

NOTE

CRAFTWORK



1. 未注尺寸公差: X±0.3; X.0±0.1; X.00分三种: (0<X≤20, ±0.005; 20<X≤150, ±0.01; X>150, ±0.015); X.000±0.005, 轴类用大尺寸, 孔类用小尺寸。  
 2. 相同形状标识: A.B.C.D.....  
 3. 未标注的螺钉深度等于螺丝直径的两倍  
 4. 未标注螺丝深度等于螺丝直径两倍。  
 5. 料直径垂直度: 0<长, 宽≤100 0.01;  
 6. 外形尺寸公差: 0.00<sup>-</sup>0.01; 孔尺寸公差: +0.01<sup>-</sup>0.00; 拖位尺寸公差: 托槽+0.01<sup>-</sup>0.03; 托0.00<sup>-</sup>0.02

UNITS

MM

NAME (INTENDED USE)

MAT'L

PP

PART NO. (INTENDED USE)

TITLE:

5ml 吸头

APPD:


DWG NO.:

5ML-D5-3

CHKD:

DR:

08-Aug-22



SCALE

1.000

SHEET

1/1

REV.

A

1

2

3

4

A4